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English

Manufacturer and European address:

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FALCON MEDICAL Instrument s I

Recommendations and notes

Before using FALCON MEDICAL instruments from class I, the surgeon is advised to take note of the following information:

General

- The use of FALCON MEDICAL instruments requires professional application in compliance with all surgical and safety regulations. No warranty or liability can be accepted for the safety and functional use of FALCON MEDICAL instruments if the instructions in the operating manual and these notes are not observed and followed. All instructions for use by FALCON MEDICAL are issued based on the best of our knowledge in accordance with the current state of medical science. However, this does not release the user from their responsibility as a medical specialist.
- The following points are particularly important to ensure the success of an operation using FALCON MEDICAL instruments:
 - Correct cleaning and disinfection
 - Correct sterilization
 - Correct maintenance and care
 - Regular inspectionsComplications that may arise from failure to observe these instructions for use and the expertise required for the operation are the responsibility of the surgeon and the specialists entrusted with the follow-up treatment and cannot be attributed to either the manufacturer or the supplier of FALCON MEDICAL products.
- The surgeon must be entirely familiar with the joint replacement procedure, the FALCON MEDICAL instruments used and the surgical technique to be followed. Written surgical instructions can be obtained from FALCON MEDICAL at any time. However, they do not claim to be exhaustive and are not sufficient for the immediate use of the tools. Instruction in the general surgical technique by an experienced surgeon is essential.
- Before use, the surgeon must ensure that the instruments are working properly.
- The user and/or patient is required to report any serious incidents that have occurred in relation to the product to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.
- FALCON MEDICAL instruments are associated parts of Falcon Medical hip implants and may only be used in conjunction with these. Additional specific instructions are contained in the operating instructions of the product.
- In the event that FALCON MEDICAL instruments are combined with third-party products, the user must take extra care to ensure the compatibility of the connections. Due to the large number of products on the market, it is impossible for FALCON MEDICAL to give a general approval of the combinability of its products.
- In the case in question and if in doubt, the user must consult the manufacturer for the purpose of ascertaining compatibility, providing details of the required product combination. FALCON MEDICAL assumes no warranty or liability for compatibility with third-party products

Symbols, intended use and materials used on the label: **see bottom of page**

Classification

All instruments offered are class I in accordance with Regulation (EU) 2017/745.

Indications, contraindications and risk factors of FALCON MEDICAL hip prostheses
(The non-invasive instruments made of suitable plastic and/or metal are used when a hip prosthesis is indicated. Therefore, the indications, contraindications and risk factors correspond to those for implanting a total hip endoprosthesis and can be found in the package insert of the corresponding hip implants.)

General

The indications, contraindications and risk factors can be relative with regard to artificial joint replacement and must always be considered in relation to the patient’s overall condition. Since all instruments are accessories for hip prostheses, the indications correspond to those of a hip prosthesis:

Indications

- The indications are:
- Advanced functional impairment of the hip joint due to degenerative, post-traumatic osteoarthritis, rheumatoid arthritis, avascular necrosis
 - Revision surgery resulting from an unsuccessful endoprosthesis (see MiniMIS stem contraindication)
 - Hip dislocation resulting from a fracture or non-reducible fractures for which adequate fixation cannot be achieved.

General contraindications

- The contraindications are:
- Acute, chronic or latent infection, locally (especially in the lower extremities) or systemic
 - Severe muscle, nerve or vascular diseases which would endanger the affected extremity
 - Inadequate bone structures which can endanger the anchoring and support of the implant
 - Any concomitant diseases which can endanger the function of the implant or which carry the risk of a life-threatening post-operative complication.
 - Known pronounced allergy to any of the material component of the implants

Additional contraindications for the MiniMIS stem

- Malalignment of the femoral neck
- Revision surgery resulting from an unsuccessful endoprosthesis with insufficient proximal anchorage

Additional contraindications for the revision stem

- Achievable anchoring distance < 10 cm due to the risk of the implant fracturing

Additional contraindications for the Monocon MIS stem

- Dorr type A femur

Risk factors

Risks which could impair the success of an operation include, in particular:

- Significant osteoporosis or osteomalacia
- Severe malformations, congenital hip dislocation
- Local bone tumours
- Systemic illnesses and metabolic disorders (especially in connection with cortisone and cytostatics)
- Anamnestic indications of infections and falls
- Patient obesity (body weight > 100 kg)
- Patient involved in intensive physical activity (e.g. running, lifting heavy weights or competitive sports)
- Allergy to any of the material components of the implant
- Muscle weakness
- Neurological diseases which could pose a functional risk for the artificial joint replacement; epilepsy, for example
- Patient who is not able to understand and follow the physician’s instructions
- Poor implantation, rehabilitation and complications.
- Patient involved in competitive sports
- Inadequate stability and joint tension of the implants (can lead to microseparation, dislocation and/or excessive stress on the tribological components leading to fracture of the components)
- Nutritional deficiencies
- Excessive nicotine consumption (smoking)
- Alcoholism or substance abuse

Possible complications
(the complications related to implanting a total hip endoprosthesis can be found in the accompanying text for the associated hip implants; only the specific complications of the instruments are listed here).

- Patient-related complications:**
- Infections
 - Extended operation time (including switching to another implant system during the operation)

Product-related complications:

- Instrument fracture (including due to excessive force)
- Damaged instrument
- Instrument stops working

Cleaning, disinfection, sterilization, maintenance, inspection and care

FALCON MEDICAL instruments have been checked for cleaning, disinfection and sterilization; the reference parameters for this can be found in the preparation instructions provided by Falcon Medical. The current version of these instructions can be obtained directly from FALCON MEDICAL. They contain sample images, reference chemicals, machine and programme parameters as well as maintenance instructions and inspection recommendations.

Foreseeable service life of the instruments

It is not possible to define a specific foreseeable service life as this depends on various factors, and the service life of the instruments depends on how well they are used, cared for and maintained. Cutting instruments (rasps, reamers) in particular can cause damage to plastic components and jeopardize their function due to improper use and carelessness. These instruments must be replaced by the user. Plastic parts and O-rings become brittle after frequent sterilization and must be replaced in accordance with the requirements set out in the preparation instructions. The following applies in this case: If cared for properly, the instruments can be reused until the end of their service life (see “Special instructions” section in RHW R6 05/2023), provided that they have not been reconditioned by the manufacturer in the meantime; any further use beyond this or the use of damaged and/or contaminated instruments is the responsibility of the user. In the event of non-compliance, any liability is excluded.

Packaging and sterility

Instruments manufactured, supplied and provided by FALCON MEDICAL are supplied in a clean but non-sterile condition either in trays or in trays and sterile containers. Individual instruments to be sent are always delivered in the non-sterile state in protective packaging, which is used exclusively for protecting them during transport. This packaging may not be used for cleaning, disinfection and sterilization.

Storage and handling of the instruments

The instruments can be stored sterile ready for use or in a non-sterile condition. The instruments may only be used in a prepared and sterilized state. The notes in the preparation instructions provided and the associated operating instructions apply.

Instructions for proper removal, disposal and return

The instruments may be owned by the customer or by Falcon Medical. FALCON MEDICAL cannot provide specific information about the proper disposal of damaged components as FALCON MEDICAL has no information about individual contamination.

Instruments owned by Falcon Medical must be prepared using a validated procedure, sterilized and returned to Falcon Medical with proof that this procedure has been completed.

Limitation of liability

Claims for damages against FALCON MEDICAL are excluded, unless they are related to compensation for personal injury, and FALCON MEDICAL is responsible for intent or gross negligence. Any claims for damages must be asserted against the manufacturer in court within six months, otherwise they shall be excluded.

Further information

Further information on the use of this product can be obtained from the manufacturer or distributor of FALCON MEDICAL hip rasps.

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