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**English**

Manufacturer and European address:

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**FALCON MEDICAL Instruments, Ila Rasps**

**Recommendations and notes**

Before using FALCON MEDICAL rasps (instruments) from class Ila, the surgeon is advised to take note of the following information:

**General**

- The use of FALCON MEDICAL instruments requires professional application in compliance with all surgical and safety regulations. No warranty or liability can be accepted for the safety and functional use of FALCON MEDICAL instruments if the instructions in the operating manual and these notes are not observed and followed. All instructions for use by FALCON MEDICAL are issued based on the best of our knowledge in accordance with the current state of medical science. However, this does not release the user from their responsibility as a medical specialist.
- The following points are particularly important to ensure the success of an operation using FALCON MEDICAL instruments:
  - Correct cleaning and disinfection
  - Correct sterilization
  - Correct maintenance and care
  - Regular inspections
Complications that may arise from failure to observe these instructions for use and the expertise required for the operation are the responsibility of the surgeon and the specialists entrusted with the follow-up treatment and cannot be attributed to either the manufacturer or the supplier of FALCON MEDICAL products.
- The surgeon must be entirely familiar with the joint replacement procedure, the FALCON MEDICAL instruments used and the surgical technique to be followed. Written surgical instructions can be obtained from FALCON MEDICAL at any time. However, they do not claim to be exhaustive and are not sufficient for the immediate use of the tools. Instruction in the general surgical technique by an experienced surgeon is essential.
- Before use, the surgeon must ensure that the instruments are working properly.
- The user and/or patient is required to report any serious incidents that have occurred in relation to the product to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.
- FALCON MEDICAL rasps are associated parts of Falcon Medical hip sockets and may only be used in conjunction with these. Additional specific instructions are contained in the operating instructions of the product.
- In the event that FALCON MEDICAL instruments are combined with third-party products, the user must take extra care to ensure the compatibility of the connections. Due to the large number of products on the market, it is impossible for FALCON MEDICAL to give a general approval of the combinability of its products.
- In the case in question and if in doubt, the user must consult the manufacturer for the purpose of ascertaining compatibility, providing details of the required product combination. FALCON MEDICAL assumes no warranty or liability for compatibility with third-party products

**Symbols and materials used on the label: see bottom of page**

**Classification**

FALCON MEDICAL rasps have been developed to support the implantation of hip prostheses. All rasps offered are class Ila in accordance with Regulation (EU) 2017/745.

**Intended use**

The intended purpose is based on the indication of the associated femoral implant. The rasp is used for internally preparing the femoral medullary canal to make sure that the corresponding implant fits securely. At the same time, an intraoperative function test can be performed with the size of rasp intended to be used.

A mechanical-pneumatic drive using a third-party product (IMT "Specht"; manufacturer IMT Integral Medizintechnik AG) is possible using the same instrument interface as for the manual preparation; this combination has been verified by Falcon Medical.

**Indications, contraindications and risk factors of FALCON MEDICAL hip prostheses**
*(Metal rasps are used when a hip prosthesis is indicat-ed. Therefore, the indications, contraindications and risk factors correspond to those for implanting a total hip endoprosthesis and can be found in the package insert of the corresponding hip implants.)*

**General**

The indications, contraindications and risk factors can be relative with regard to artificial joint replacement and must always be considered in relation to the patient's overall condition. Since all instruments are accessories for hip prostheses, the indications correspond to those of a hip prosthesis:

**Indications**

The indications are:

- a) Advanced functional impairment of the hip joint due to degenerative, post-traumatic osteoarthritis, rheumatoid arthritis, avascular necrosis
- b) Revision surgery resulting from an unsuccessful endoprosthesis (see MinimIS stem contraindication)
- c) Hip dislocation resulting from a fracture or non-reducible fractures for which adequate fixation cannot be achieved.

**General contraindications**

- a) Acute, chronic or latent infection, locally (especially in the lower extremities) or systemic
- b) Severe muscle, nerve or vascular diseases which would endanger the affected extremity
- c) Inadequate bone structures which can endanger the anchoring and support of the implant
- d) Any concomitant diseases which can endanger the function of the implant or which carry the risk of a life-threatening post-operative complication.
- e) Known pronounced allergy to any of the material component of the implants

**Additional contraindications for the MinimIS stem**

- a) Malalignment of the femoral neck
- b) Revision surgery resulting from an unsuccessful endoprosthesis with insufficient proximal anchorage

**Additional contraindications for the revision stem**

- a) Achievable anchoring distance < 10 cm due to the risk of the implant fracturing

**Additional contraindications for the Monoco MIS stem**

- a) Dorr type A femur

**Risk factors**

Risks which could impair the success of an operation include, in particular:

- a) Significant osteoporosis or osteomalacia
- b) Severe malformations, congenital hip dislocation
- c) Local bone tumours
- d) Systemic illnesses and metabolic disorders (especially in connection with cortisone and cytostatics)
- e) Anaesthetic indications of infections and falls
- f) Patient obesity (body weight > 100 kg)
- g) Patient involved in intensive physical activity (e.g. running, lifting heavy weights or competitive sports)
- h) Allergy to any of the material components of the implant
- i) Muscle weakness
- j) Neurological diseases which could pose a functional risk for the artificial joint replacement; epilepsy, for example
- k) Patient who is not able to understand and follow the physician's instructions
- l) Poor implantation, rehabilitation and complications.
- m) Patient involved in competitive sports
- n) Inadequate stability and joint tension of the implants (can lead to microseparation, dislocation and excessive stress on the tribological components leading to fracture of the components)
- o) Nutritional deficiencies
- p) Excessive nicotine consumption (smoking)
- q) Alcoholism or substance abuse

**Possible complications**
*(the complications related to implanting a total hip endoprosthesis can be found in the accompanying text for the associated hip implants; only the specific complications of the rasps are listed here.)*

**Patient-related complications:**

- a) Bone injuries (chipping or even fractures)
- b) Soft tissue injuries (including nerves, vessels and muscles)
- c) Infections
- d) Extended operation time (including switching to another implant system during the operation)
- e) Tissue reaction (allergy, sensitization, iritation, etc.)

**Product-related complications:**

- a) Instrument fracture (including due to excessive force)
- b) Damaged instrument
- c) Instrument stops working
- d) Holder becomes mishapen (e.g. rasp adapter, threaded cup holder)

**Cleaning, disinfection, sterilization, maintenance, inspection and care**

FALCON MEDICAL instruments have been checked for cleaning, disinfection and sterilization; the refer-ence parameters for this can be found in the preparation instructions provided by Falcon Medical. The current version of these instructions can be obtained directly from FALCON MEDICAL. They contain sample images, reference chemicals, machine and programme parameters as well as maintenance instructions and inspection recommendations.

**Foreseeable service life of the instruments**

It is not possible to define a specific foreseeable service life as this depends on various factors, and the service life of the instruments depends on how well they are used, cared for and maintained. Special cutting instruments (rasps) in particular can be damaged through improper use and carelessness, and this can impair their function. These instruments must be replaced by the user. Plastic parts and O-rings become brittle after frequent sterilization and must be replaced in accordance with the requirements set out in the preparation instructions.

The following applies in this case:

If cared for properly, the instruments can be reused until the end of their service life (see "Special instru-ctions" section in RHW R6 05/2023), provided that they have not been reconditioned by the manufacturer in the meantime; any further use beyond this or the use of damaged and/or contaminated instruments is the responsibility of the user.

In the event of non-compliance, any liability is excluded.

**Packaging and sterility**

Instrument manufactured, supplied and provided by FALCON MEDICAL are supplied in a clean but non-sterile condition either in trays or in trays and sterile containers. Individual instruments to be sent are always delivered in the non-sterile state in protective packaging, which is used exclusively for protecting them during transport. This packaging may not be used for cleaning, disinfection and sterilization.

**Storage and handling of the instruments**

The instruments can be stored sterile ready for use or in a non-sterile condition. The instruments may only be used in a prepared and sterilized state. The notes in the preparation instructions provided and the associated operating instructions apply.

**Instructions for proper removal, disposal and return**

The instruments may be owned by the customer or by Falcon Medical. FALCON MEDICAL cannot pro-vide specific information about the proper disposal of damaged components as FALCON MEDICAL has no information about individual contamination.

Instruments owned by Falcon Medical must be prepared using a validated procedure, sterilized and returned to Falcon Medical with proof that this procedure has been completed.

**Limitation of liability**

Claims for damages against FALCON MEDICAL are excluded, unless they are related to compensation for personal injury, or FALCON MEDICAL is responsible for intent or gross negligence. Any claims for damages must be asserted against the manufacturer in court within six months, otherwise they shall be excluded.

**Further information**

Further information on the use of this product can be obtained from the manufacturer or distributor of FAL-CON MEDICAL hip rasps.

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